

## LYPLAL1

---

**Reactivity:**Human Mouse

**Tested applications:**WB IHC

**Recommended Dilution:**WB 1:200 - 1:2000 IHC 1:20 - 1:200

**Calculated MW:**26kDa

**Observed MW:**Refer to figures

**Immunogen:**

A synthetic Peptide of human LYPLAL1

**Storage Buffer:**

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

Q96AV0;

*To place an order, please [Click HERE](#).*

**Catalog #:**A5135

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**127018

**Isotype:**IgG

**Swiss Prot:**Q5VWZ2

**Purity:**Affinity purification

For research use only.