

## COMMD1

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB IHC IP

**Recommended Dilution:**WB 1:500 - 1:2000 IHC 1:20 - 1:200 IP 1:20 - 1:50

**Calculated MW:**21kDa

**Observed MW:**Refer to figures

**Immunogen:**

A synthetic Peptide of human COMMD1

**Storage Buffer:**

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

MURR1; C2orf5;

**Catalog #:**A5174

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**150684

**Isotype:**IgG

**Swiss Prot:**Q8N668

**Purity:**Affinity purification

For research use only.

**Background:**

COMMD1 is a regulator of copper homeostasis, sodium uptake, and NF-kappa-B (see MIM 164011) signaling (de Bie et al., 2005

*To place an order, please [Click HERE](#).*