

JAKMIP1

Reactivity: Human

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 73kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human JAKMIP1

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

JAMIP1; MARLIN1; Gababrbp;

To place an order, please [Click HERE](#).

Catalog #: A5178

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 152789

Isotype: IgG

Swiss Prot: Q96N16

Purity: Affinity purification

For research use only.