

RASSF6

Reactivity:Mouse Rat

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:43kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human RASSF6

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

n

Catalog #:A5187

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:166824

Isotype:IgG

Swiss Prot:Q6ZTQ3

Purity:Affinity purification

For research use only.

Background:

This gene encodes a member of the Ras-association domain family (RASSF). Members of this family form the core of a highly conserved tumor suppressor network, the Salvador-Warts-Hippo (SWH) pathway. The protein encoded by this gene is a Ras effector protein that induces apoptosis. A genomic region containing this gene has been linked to susceptibility to viral bronchiolitis. Alternative splicing results in multiple transcript variants and protein isoforms.

To place an order, please [Click HERE](#).