

RASSF6

Reactivity: Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 43kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human RASSF6

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

n

Catalog #: A5187

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 166824

Isotype: IgG

Swiss Prot: Q6ZTQ3

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the Ras-association domain family (RASSF). Members of this family form the core of a highly conserved tumor suppressor network, the Salvador-Warts-Hippo (SWH) pathway. The protein encoded by this gene is a Ras effector protein that induces apoptosis. A genomic region containing this gene has been linked to susceptibility to viral bronchiolitis. Alternative splicing results in multiple transcript variants and protein isoforms.

To place an order, please [Click HERE](#).