

USP12

Reactivity:Human Mouse Rat

Tested applications:WB

Recommended Dilution:WB 1:200 - 1:2000

Calculated MW:43kDa

Observed MW:Refer to figures

Immunogen:

A synthetic Peptide of human USP12

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

UBH1; USP12L1;

To place an order, please [Click HERE](#).

Catalog #:A5201

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:219333

Isotype:IgG

Swiss Prot:O75317

Purity:Affinity purification

For research use only.