

## TMED4

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB IF

**Recommended Dilution:**WB 1:500 - 1:2000 IF 1:10 - 1:100

**Calculated MW:**26kDa

**Observed MW:**Refer to figures

**Immunogen:**

A synthetic Peptide of human TMED4

**Storage Buffer:**

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

HNLF; ERS25; GMP25iso; p24alpha3;

*To place an order, please [Click HERE](#).*

**Catalog #:**A5208

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**222068

**Isotype:**IgG

**Swiss Prot:**Q7Z7H5

**Purity:**Affinity purification

For research use only.