

## LCORL

---

**Reactivity:** Human Mouse Rat

**Tested applications:** WB IHC

**Recommended Dilution:** WB 1:500 - 1:2000 IHC 1:20 - 1:200

**Calculated MW:** 67kDa

**Observed MW:** Refer to figures

**Immunogen:**

A synthetic Peptide of human LCORL

**Storage Buffer:**

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

MLR1;

**Catalog #:** A5209

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 254251

**Isotype:** IgG

**Swiss Prot:** Q8N3X6

**Purity:** Affinity purification

For research use only.

**Background:**

This gene encodes a transcription factor that appears to function in spermatogenesis.

Polymorphisms in this gene are associated with measures of skeletal frame size and adult height.

Alternative splicing results in multiple transcript variants.

*To place an order, please [Click HERE](#).*