

IAH1

Reactivity:Human Mouse Rat

Tested applications:WB IP

Recommended Dilution:WB 1:500 - 1:2000 IP 1:20 - 1:50

Calculated MW:28kDa

Observed MW:Refer to figures

Immunogen:

A synthetic Peptide of human IAH1

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

To place an order, please [Click HERE](#).

Catalog #:A5217

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:285148

Isotype:IgG

Swiss Prot:Q2TAA2

Purity:Affinity purification

For research use only.