

MMAB

Reactivity:Human Mouse

Tested applications:WB IHC IP

Recommended Dilution:WB 1:200 - 1:2000 IHC 1:20 - 1:200 IP 1:20 - 1:50

Calculated MW:27kDa

Observed MW:Refer to figures

Immunogen:

A synthetic Peptide of human MMAB

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ATR; cob; cbIB; CFAP23;

Catalog #:A5221

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:326625

Isotype:IgG

Swiss Prot:Q96EY8

Purity:Affinity purification

For research use only.

Background:

This gene encodes a protein that catalyzes the final step in the conversion of vitamin B(12) into adenosylcobalamin (AdoCbl), a vitamin B12-containing coenzyme for methylmalonyl-CoA mutase. Mutations in the gene are the cause of vitamin B12-dependent methylmalonic aciduria linked to the cbIB complementation group. Alternatively spliced transcript variants have been found.

To place an order, please [Click HERE](#).