

## CISD2

---

**Reactivity:**Human Rat

**Tested applications:**WB

**Recommended Dilution:**WB 1:500 - 1:2000

**Calculated MW:**15kDa

**Observed MW:**Refer to Figures

**Immunogen:**

Recombinant protein of human CISD2

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Concentration:**

j

**Synonym:**

ERIS; WFS2; ZCD2; NAF-1; Miner1;

**Catalog #:**A5231

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**493856

**Isotype:**IgG

**Swiss Prot:**Q8N5K1

**Purity:**Affinity purification

For research use only.

**Background:**

The protein encoded by this gene is a zinc finger protein that localizes to the endoplasmic reticulum. The encoded protein binds an iron/sulfur cluster and may be involved in calcium homeostasis. Defects in this gene are a cause of Wolfram syndrome 2.

*To place an order, please [Click HERE](#).*