

CPLX3

Reactivity:Human Mouse Rat

Tested applications:WB IF IP

Recommended Dilution:WB 1:200 - 1:2000 IF 1:10 - 1:100 IP 1:20 - 1:50

Calculated MW:18kDa

Observed MW:Refer to figures

Immunogen:

A synthetic Peptide of human CPLX3

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CPXIII; CPX-III; Nbla11589;

To place an order, please [Click HERE](#).

Catalog #:A5234

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:594855

Isotype:IgG

Swiss Prot:Q8WVH0

Purity:Affinity purification

For research use only.