

## GTF2H2C

---

**Reactivity:**Human

**Tested applications:**WB

**Recommended Dilution:**WB 1:200 - 1:2000

**Calculated MW:**44kDa

**Observed MW:**Refer to figures

**Immunogen:**

Recombinant protein of human GTF2H2C

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Concentration:**

e

**Synonym:**

GTF2H2C\_2;

*To place an order, please [Click HERE](#).*

**Catalog #:**A5239

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**728340

**Isotype:**IgG

**Swiss Prot:**Q6P1K8

**Purity:**Affinity purification

For research use only.