

## NRF1

---

**Reactivity:**Human

**Tested applications:**WB

**Recommended Dilution:**WB 1:500 - 1:2000

**Calculated MW:**54kDa

**Observed MW:**Refer to Figures

**Immunogen:**

Recombinant protein of human NRF1

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Concentration:**

hi

**Synonym:**

NRF1; ALPHA-PAL; Nuclear respiratory factor 1; Alpha palindromic-binding protein;

**Catalog #:**A5264

**Antibody Type:**

Monoclonal Antibody

**Species:**Mouse

**Gene ID:**4899

**Isotype:**IgG

**Swiss Prot:**Q16656

**Purity:**Affinity purification

For research use only.

**Background:**

This gene encodes a protein that homodimerizes and functions as a transcription factor which activates the expression of some key metabolic genes regulating cellular growth and nuclear genes required for respiration, heme biosynthesis, and mitochondrial DNA transcription and replication. The protein has also been associated with the regulation of neurite outgrowth. Alternate transcriptional splice variants, which encode the same protein, have been characterized. Additional variants encoding different protein isoforms have been described but they have not been fully characterized. Confusion has occurred in bibliographic databases due to the shared symbol of NRF1 for this gene and for "nuclear factor (erythroid-derived 2)-like 1" which has an official symbol of NFE2L1.

**To place an order, please [Click HERE](#).**