

EIF4A1

Reactivity:Human Mouse Rat

Tested applications:WB IF

Recommended Dilution:WB 1:500 - 1:2000 IF 1:20 - 1:50

Calculated MW:46kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human EIF4A1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

j

Synonym:

DDX2A; EIF4A; EIF-4A; eIF4A-I; eIF-4A-I

Catalog #:A5294

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:1973

Isotype:IgG

Swiss Prot:P60842

Purity:Affinity purification

For research use only.

Background:

EIF4A1 is an ATP-dependent RNA helicase which is a subunit of the eIF4F complex involved in cap recognition and is required for mRNA binding to ribosome. In the current model of translation initiation, eIF4A unwinds RNA secondary structures in the 5'-UTR of mRNAs which is necessary to allow efficient binding of the small ribosomal subunit, and subsequent scanning for the initiator codon.

To place an order, please [Click HERE](#).