

## AIMP1

**Reactivity:** Human Mouse Rat

**Tested applications:** WB IHC IF

**Recommended Dilution:** WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

**Calculated MW:** 34kDa

**Observed MW:** Refer to Figures

**Immunogen:**

Recombinant protein of human AIMP1

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

p43; HLD3; EMAP2; SCYE1; EMAPII;

**Catalog #:** A5301

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 9255

**Isotype:** IgG

**Swiss Prot:** Q12904

**Purity:** Affinity purification

For research use only.

**Background:**

The protein encoded by this gene is a cytokine that is specifically induced by apoptosis, and it is involved in the control of angiogenesis, inflammation, and wound healing. The release of this cytokine renders the tumor-associated vasculature sensitive to tumor necrosis factor. The precursor protein is identical to the p43 subunit, which is associated with the multi-tRNA synthetase complex, and it modulates aminoacylation activity of tRNA synthetase in normal cells. This protein is also involved in the stimulation of inflammatory responses after proteolytic cleavage in tumor cells. Multiple transcript variants encoding different isoforms have been found for this gene. A pseudogene has been identified on chromosome 20.

**To place an order, please [Click HERE](#).**