

## NDRG2

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB IHC

**Recommended Dilution:**WB 1:500 - 1:2000 IHC 1:50 - 1:200

**Calculated MW:**41kDa

**Observed MW:**Refer to Figures

**Immunogen:**

Recombinant protein of human NDRG2

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

SYLD;

**Catalog #:**A5319

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**57447

**Isotype:**IgG

**Swiss Prot:**Q9UN36

**Purity:**Affinity purification

For research use only.

**Background:**

This gene is a member of the N-myc downregulated gene family which belongs to the alpha/beta hydrolase superfamily. The protein encoded by this gene is a cytoplasmic protein that may play a role in neurite outgrowth. This gene may be involved in glioblastoma carcinogenesis. Several alternatively spliced transcript variants of this gene have been described, but the full-length nature of some of these variants has not been determined.

*To place an order, please [Click HERE](#).*