

NDRG2

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:41kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human NDRG2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

SYLD;

Catalog #:A5319

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:57447

Isotype:IgG

Swiss Prot:Q9UN36

Purity:Affinity purification

For research use only.

Background:

This gene is a member of the N-myc downregulated gene family which belongs to the alpha/beta hydrolase superfamily. The protein encoded by this gene is a cytoplasmic protein that may play a role in neurite outgrowth. This gene may be involved in glioblastoma carcinogenesis. Several alternatively spliced transcript variants of this gene have been described, but the full-length nature of some of these variants has not been determined.

To place an order, please [Click HERE](#).