

NOV

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW: 39kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human NOV

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CCN3; NOVh; IBP-9; IGFBP9; IGFBP-9;

Catalog #: A5320

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 4856

Isotype: IgG

Swiss Prot: P48745

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a small secreted cysteine-rich protein and a member of the CCN family of regulatory proteins. CNN family proteins associate with the extracellular matrix and play an important role in cardiovascular and skeletal development, fibrosis and cancer development.

To place an order, please [Click HERE](#).