

## NOV

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB IHC

**Recommended Dilution:**WB 1:500 - 1:2000 IHC 1:50 - 1:100

**Calculated MW:**39kDa

**Observed MW:**Refer to figures

**Immunogen:**

Recombinant protein of human NOV

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

CCN3; NOVh; IBP-9; IGFBP9; IGFBP-9;

**Background:**

The protein encoded by this gene is a small secreted cysteine-rich protein and a member of the CCN family of regulatory proteins. CNN family proteins associate with the extracellular matrix and play an important role in cardiovascular and skeletal development, fibrosis and cancer development.

*To place an order, please [Click HERE](#).*

**Catalog #:**A5320

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**4856

**Isotype:**IgG

**Swiss Prot:**P48745

**Purity:**Affinity purification

For research use only.