

PTGES3

Reactivity:Human Mouse Rat

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:20 - 1:50

Calculated MW:19kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human PTGES3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

P23; TEBP; cPGES;

To place an order, please [Click HERE](#).

Catalog #:A5325

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:10728

Isotype:IgG

Swiss Prot:Q15185

Purity:Affinity purification

For research use only.