

BPI

Reactivity:Human Mouse

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW:54kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human BPI

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

rBPI; BPIFD1

Catalog #:A5338

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:671

Isotype:IgG

Swiss Prot:P17213

Purity:Affinity purification

For research use only.

Background:

This gene encodes a lipopolysaccharide binding protein. It is associated with human neutrophil granules and has bactericidal activity on gram-negative organisms.

To place an order, please [Click HERE](#).