

HDGF

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF IP RIP

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200 IP 1:20 - 1:50
RIP 1:20 - 1:50

Calculated MW: 26kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human HDGF

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

j

Synonym:

HMG1L2;

Catalog #: A5347

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 3068

Isotype: IgG

Swiss Prot: P51858

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the hepatoma-derived growth factor family. The encoded protein has mitogenic and DNA-binding activity and may play a role in cellular proliferation and differentiation. This gene was thought initially to be located on chromosome X, however, that location has been determined to correspond to a related pseudogene. Alternatively spliced transcript variants encoding distinct isoforms have been described.

To place an order, please [Click HERE](#).