

VAPB

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:20 - 1:50

Calculated MW: 27kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human VAPB

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ALS8; VAP-B; VAMP-B;

Catalog #: A5363

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 9217

Isotype: IgG

Swiss Prot: O95292

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a type IV membrane protein found in plasma and intracellular vesicle membranes. The encoded protein is found as a homodimer and as a heterodimer with VAPA. This protein also can interact with VAMP1 and VAMP2 and may be involved in vesicle trafficking.

To place an order, please [Click HERE](#).