

CALD1

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW: 63kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human CALD1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CDM; HCAD; LCAD; H-CAD; L-CAD; NAG22

Catalog #: A5366

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 800

Isotype: IgG

Swiss Prot: Q05682

Purity: Affinity purification

For research use only.

Background:

This gene encodes a calmodulin- and actin-binding protein that plays an essential role in the regulation of smooth muscle and nonmuscle contraction. The conserved domain of this protein possesses the binding activities to Ca(2+)-calmodulin, actin, tropomyosin, myosin, and phospholipids. This protein is a potent inhibitor of the actin-tropomyosin activated myosin MgATPase, and serves as a mediating factor for Ca(2+)-dependent inhibition of smooth muscle contraction. Alternative splicing of this gene results in multiple transcript variants encoding distinct isoforms.

To place an order, please [Click HERE](#).