

RAB3A

Reactivity:Human Mouse Rat

Tested applications:WB IHC IP

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IP 1:20 - 1:50

Calculated MW:25kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human RAB3A

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

To place an order, please [Click HERE](#).

Catalog #:A5379

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:5864

Isotype:IgG

Swiss Prot:P20336

Purity:Affinity purification

For research use only.