

SFRP2

Reactivity:Human Mouse Rat

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW:33kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human SFRP2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

FRP-2; SARP1; SDF-5;

Catalog #:A5383

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:6423

Isotype:IgG

Swiss Prot:Q96HF1

Purity:Affinity purification

For research use only.

Background:

This gene encodes a member of the SFRP family that contains a cysteine-rich domain homologous to the putative Wnt-binding site of Frizzled proteins. SFRPs act as soluble modulators of Wnt signaling. Methylation of this gene is a potential marker for the presence of colorectal cancer.

To place an order, please [Click HERE](#).