

ASH2L

Reactivity:Human Mouse Rat

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:20 - 1:50

Calculated MW:69kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human ASH2L

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

f

Synonym:

ASH2; Bre2; ASH2L1; ASH2L2;

To place an order, please [Click HERE](#).

Catalog #:A5392

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:9070

Isotype:IgG

Swiss Prot:Q9UBL3

Purity:Affinity purification

For research use only.