

C7

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 94kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human C7

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #: A5394

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 730

Isotype: IgG

Swiss Prot: P10643

Purity: Affinity purification

For research use only.

Background:

C7 is a component of the complement system. It participates in the formation of Membrane Attack Complex (MAC). People with C7 deficiency are prone to bacterial infection.

To place an order, please [Click HERE](#).