

CARTPT

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 13kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human CARTPT

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CART;

Catalog #: A5395

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 9607

Isotype: IgG

Swiss Prot: Q16568

Purity: Affinity purification

For research use only.

Background:

This gene encodes a secreted protein which is processed by prohormone/proprotein convertases to produce smaller, biologically active peptides. Expression of the transcript for this gene is regulated by certain drugs such as cocaine, and the encoded protein is thought to be involved in the regulation of appetite and stress. Mutations in this gene are associated with susceptibility to obesity.

To place an order, please [Click HERE](#).