

MPG

Reactivity:Human Mouse Rat

Tested applications:WB IHC IF IP

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100 IP 1:50 - 1:100

Calculated MW:32kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human MPG

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

ef

Synonym:

AAG; MDG; ADPG; APNG; Mid1; anpg; PIG11; PIG16; CRA36.1;

To place an order, please [Click HERE](#).

Catalog #:A5410

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:4350

Isotype:IgG

Swiss Prot:P29372

Purity:Affinity purification

For research use only.