

NPC2

Reactivity:Human Mouse Rat

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW:17kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human NPC2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

HE1; EDDM1;

Background:

This gene encodes a protein containing a lipid recognition domain. The encoded protein may function in regulating the transport of cholesterol through the late endosomal/lysosomal system. Mutations in this gene have been associated with Niemann-Pick disease, type C2 and frontal lobe atrophy.

To place an order, please [Click HERE](#).

Catalog #:A5413

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:10577

Isotype:IgG

Swiss Prot:P61916

Purity:Affinity purification

For research use only.