

## TFDP1

**Reactivity:** Human Mouse Rat

**Tested applications:** WB IHC

**Recommended Dilution:** WB 1:500 - 1:2000 IHC 1:50 - 1:200

**Calculated MW:** 45kDa

**Observed MW:** Refer to Figures

**Immunogen:**

Recombinant protein of human TFDP1

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Concentration:**

bk

**Synonym:**

DP1; Dp-1; DRTF1;

**Catalog #:** A5422

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 7027

**Isotype:** IgG

**Swiss Prot:** Q14186

**Purity:** Affinity purification

For research use only.

**Background:**

This gene encodes a member of a family of transcription factors that heterodimerize with E2F proteins to enhance their DNA-binding activity and promote transcription from E2F target genes. The encoded protein functions as part of this complex to control the transcriptional activity of numerous genes involved in cell cycle progression from G1 to S phase. Alternative splicing results in multiple transcript variants. Pseudogenes of this gene are found on chromosomes 1, 15, and X.

*To place an order, please [Click HERE](#).*