

## Suz12

---

**Reactivity:**Human Mouse

**Tested applications:**WB

**Recommended Dilution:**WB 1:500 - 1:2000

**Calculated MW:**83kDa

**Observed MW:**Refer to Figures

**Immunogen:**

A synthetic peptide of human Suz12

**Storage Buffer:**

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Concentration:**

f

**Synonym:**

CHET9; JJAZ1;

**Catalog #:**A5429

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**23512

**Isotype:**IgG

**Swiss Prot:**Q15022

**Purity:**Affinity purification

For research use only.

**Background:**

This zinc finger gene has been identified at the breakpoints of a recurrent chromosomal translocation reported in endometrial stromal sarcoma. Recombination of these breakpoints results in the fusion of this gene and JAZF1. The protein encoded by this gene contains a zinc finger domain in the C terminus of the coding region.

*To place an order, please [Click HERE](#).*