

## TMSB4X

---

**Reactivity:** Human Mouse Rat

**Tested applications:** WB IHC IF

**Recommended Dilution:** WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

**Calculated MW:** 5kDa

**Observed MW:** Refer to figures

**Immunogen:**

Recombinant protein of human TMSB4X

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

FX; TB4X; PTMB4; TMSB4;

**Catalog #:** A5438

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 7114

**Isotype:** IgG

**Swiss Prot:** P62328

**Purity:** Affinity purification

For research use only.

**Background:**

This gene encodes an actin sequestering protein which plays a role in regulation of actin polymerization. The protein is also involved in cell proliferation, migration, and differentiation. This gene escapes X inactivation and has a homolog on chromosome Y.

*To place an order, please [Click HERE](#).*