

PI3

Reactivity:Human Mouse

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:100 - 1:200

Calculated MW:12kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human PI3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

Cementoin; Elafin; Elastase specific inhibitor; ESI; Peptidase inhibitor 3; PI 3; PI3; Protease inhibitor WAP3; SKALP; WAP3; WFDC14

To place an order, please [Click HERE](#).

Catalog #:A5441

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:5266

Isotype:IgG

Swiss Prot:P19957

Purity:Affinity purification

For research use only.