

CBR1

Reactivity:Human Mouse Rat

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:10 - 1:100

Calculated MW:30kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human CBR1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CBR; hCBR1; SDR21C1

Catalog #:A5446

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:873

Isotype:IgG

Swiss Prot:P16152

Purity:Affinity purification

For research use only.

Background:

The protein encoded by this gene belongs to the short-chain dehydrogenases/reductases (SDR) family, which function as NADPH-dependent oxidoreductases having wide specificity for carbonyl compounds, such as quinones, prostaglandins, and various xenobiotics. Alternatively spliced transcript variants have been found for this gene.

To place an order, please [Click HERE](#).