

EIF3E

Reactivity:Human Mouse Rat

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW:52kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human EIF3E

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

INT6; EIF3S6; EIF3-P48; eIF3-p46;

To place an order, please [Click HERE](#).

Catalog #:A5447

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:3646

Isotype:IgG

Swiss Prot:P60228

Purity:Affinity purification

For research use only.