

CDC34

Reactivity:Human Mouse Rat

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW:34kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human CDC34

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

fl

Synonym:

UBC3; UBCH3; UBE2R1; E2-CDC34

Catalog #:A5457

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:997

Isotype:IgG

Swiss Prot:P49427

Purity:Affinity purification

For research use only.

Background:

The protein encoded by this gene is a member of the ubiquitin-conjugating enzyme family.

Ubiquitin-conjugating enzyme catalyzes the covalent attachment of ubiquitin to other proteins. This protein is a part of the large multiprotein complex, which is required for ubiquitin-mediated degradation of cell cycle G1 regulators, and for the initiation of DNA replication.

To place an order, please [Click HERE](#).