

MB

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW: 17kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human MB

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

PVALB; myoglobin;

Catalog #: A5471

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 4151

Isotype: IgG

Swiss Prot: P02144

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the globin superfamily and is expressed in skeletal and cardiac muscles. The encoded protein is a haemoprotein contributing to intracellular oxygen storage and transcellular facilitated diffusion of oxygen. At least three alternatively spliced transcript variants encoding the same protein have been reported.

To place an order, please [Click HERE](#).