

PRPF3

Reactivity:Human Mouse Rat

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW:78kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human PRPF3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

PRP3; RP18; HPRP3; Prp3p; HPRP3P; SNRNP90;

Background:

The removal of introns from nuclear pre-mRNAs occurs on complexes called spliceosomes, which are made up of 4 small nuclear ribonucleoprotein (snRNP) particles and an undefined number of transiently associated splicing factors. This gene product is one of several proteins that associate with U4 and U6 snRNPs. Mutations in this gene are associated with retinitis pigmentosa-18.

To place an order, please [Click HERE](#).

Catalog #:A5482

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:9129

Isotype:IgG

Swiss Prot:O43395

Purity:Affinity purification

For research use only.