

DNAJB1

Reactivity:Human Mouse Rat

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW:38kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human DNAJB1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

Hdj1; Sis1; HSPF1; Hsp40; RSPH16B;

To place an order, please [Click HERE](#).

Catalog #:A5504

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:3337

Isotype:IgG

Swiss Prot:P25685

Purity:Affinity purification

For research use only.