

## DLL1

---

**Reactivity:**Human

**Tested applications:**WB

**Recommended Dilution:**WB 1:500 - 1:2000

**Calculated MW:**78kDa

**Observed MW:**Refer to Figures

**Immunogen:**

Recombinant protein of human DLL1

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Concentration:**

t

**Synonym:**

DL1; Delta; DELTA1;

**Catalog #:**A5509

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**28514

**Isotype:**IgG

**Swiss Prot:**O00548

**Purity:**Affinity purification

For research use only.

**Background:**

DLL1 is a human homolog of the Notch Delta ligand and is a member of the delta/serrate/jagged family. It plays a role in mediating cell fate decisions during hematopoiesis. It may play a role in cell-to-cell communication.

*To place an order, please [Click HERE](#).*