

TWIST1

Reactivity:Human Mouse Rat

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:1000

Calculated MW:20kDa

Observed MW:Refer to Figures

Immunogen:

A synthetic peptide of human TWIST1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

fk

Synonym:

SCS;ACS3;CRS1;BPES2;BPES3;

Catalog #:A5526

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:7291

Isotype:IgG

Swiss Prot:Q15672

Purity:Affinity purification

For research use only.

Background:

Basic helix-loop-helix (bHLH) transcription factors have been implicated in cell lineage determination and differentiation. The protein encoded by this gene is a bHLH transcription factor and shares similarity with another bHLH transcription factor, Dermo1. The strongest expression of this mRNA is in placental tissue; in adults, mesodermally derived tissues express this mRNA preferentially. Mutations in this gene have been found in patients with Saethre-Chotzen syndrome.

To place an order, please [Click HERE](#).