

## TWIST1

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB

**Recommended Dilution:**WB 1:500 - 1:1000

**Calculated MW:**20kDa

**Observed MW:**Refer to Figures

**Immunogen:**

A synthetic peptide of human TWIST1

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Concentration:**

fk

**Synonym:**

SCS;ACS3;CRS1;BPES2;BPES3;

**Catalog #:**A5526

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**7291

**Isotype:**IgG

**Swiss Prot:**Q15672

**Purity:**Affinity purification

For research use only.

**Background:**

Basic helix-loop-helix (bHLH) transcription factors have been implicated in cell lineage determination and differentiation. The protein encoded by this gene is a bHLH transcription factor and shares similarity with another bHLH transcription factor, Dermo1. The strongest expression of this mRNA is in placental tissue; in adults, mesodermally derived tissues express this mRNA preferentially. Mutations in this gene have been found in patients with Saethre-Chotzen syndrome.

*To place an order, please [Click HERE](#).*