

CBX6

Reactivity:Human Mouse Rat

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:44kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human CBX6

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

f

Catalog #:A5533

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:23466

Isotype:IgG

Swiss Prot:O95503

Purity:Affinity purification

For research use only.

Background:

The Chromobox domain (Cbx) gene family, consisting of Polycomb and Heterochromatin Protein 1 genes, involves in transcriptional repression, cell cycle regulation and chromatin remodeling [PMID: 17486619]. Cbx6, other named neuronal pentraxin receptor (NPR), first described as a putative integral membrane pentraxin that interacts with neuronal pentraxin 1 and 2, has also been classified as a Pc protein [PMID: 14647293]. NPR was determined to be a member of a distinct pentraxin protein subfamily, and found to form heteropentamers with NP1 and NP2, followed by release from cell membranes [PMID: 10748068].

To place an order, please [Click HERE](#).