

EZH2

Reactivity:Human Mouse

Tested applications:WB IHC IF FC

Recommended Dilution:WB 1:500 - 1:1000 IHC 1:20 - 1:50 IF 1:20 - 1:50 FC 1:20 - 1:50

Calculated MW:85kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human EZH2

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

f

Synonym:

WVS; ENX1; EZH1; KMT6; WVS2; ENX-1; EZH2b; KMT6A;

Catalog #:A5535

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:2146

Isotype:IgG

Swiss Prot:Q15910

Purity:Affinity purification

For research use only.

Background:

This gene encodes a member of the Polycomb-group (PcG) family. PcG family members form multimeric protein complexes, which are involved in maintaining the transcriptional repressive state of genes over successive cell generations. This protein associates with the embryonic ectoderm development protein, the VAV1 oncoprotein, and the X-linked nuclear protein. This protein may play a role in the hematopoietic and central nervous systems. Multiple alternatively spliced transcript variants encoding distinct isoforms have been identified for this gene.

To place an order, please [Click HERE](#).