

## SKIV2L2

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB IHC

**Recommended Dilution:**WB 1:500 - 1:2000 IHC 1:50 - 1:200

**Calculated MW:**118kDa

**Observed MW:**Refer to Figures

**Immunogen:**

Recombinant protein of human SKIV2L2

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Concentration:**

j

**Synonym:**

Dob1; Mtr4; fSAP118; KIAA0052

*To place an order, please [Click HERE](#).*

**Catalog #:**A5576

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**23517

**Isotype:**IgG

**Swiss Prot:**P42285

**Purity:**Affinity purification

For research use only.