

## DHX9

---

**Reactivity:**Human

**Tested applications:**WB IHC IF

**Recommended Dilution:**WB 1:1000 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

**Calculated MW:**141kDa

**Observed MW:**Refer to Figures

**Immunogen:**

A synthetic peptide of human DHX9

**Storage Buffer:**

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Concentration:**

j

**Synonym:**

LKP; RHA; DDX9; NDH2; NDHII;

**Catalog #:**A5585

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**1660

**Isotype:**IgG

**Swiss Prot:**Q08211

**Purity:**Affinity purification

For research use only.

**Background:**

This gene encodes a member of the DEAH-containing family of RNA helicases. The encoded protein is an enzyme that catalyzes the ATP-dependent unwinding of double-stranded RNA and DNA-RNA complexes. This protein localizes to both the nucleus and the cytoplasm and functions as a transcriptional regulator. This protein may also be involved in the expression and nuclear export of retroviral RNAs. Alternate splicing results in multiple transcript variants. Pseudogenes of this gene are found on chromosomes 11 and 13.

**To place an order, please [Click HERE](#).**