

SP1

Reactivity:Human

Tested applications:WB IHC ICC

Recommended Dilution:WB 1:1000 - 1:2000 IHC 1:100 - 1:200 ICC 1:100 - 1:200

Calculated MW:90kDa

Observed MW:Refer to Figures

Immunogen:

A synthetic peptide of human SP1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

hik

Synonym:

SP1; Transcription factor Sp1;

Catalog #:A5591

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:6667

Isotype:IgG

Swiss Prot:P08047

Purity:Affinity purification

For research use only.

Background:

The protein encoded by this gene is a zinc finger transcription factor that binds to GC-rich motifs of many promoters. The encoded protein is involved in many cellular processes, including cell differentiation, cell growth, apoptosis, immune responses, response to DNA damage, and chromatin remodeling. Post-translational modifications such as phosphorylation, acetylation, glycosylation, and proteolytic processing significantly affect the activity of this protein, which can be an activator or a repressor. Three transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).